



Parks Medical Corporation

Affordable and Innovative Endoscopy Technology

Introducing the Canady ESU and Plasma Coagulation (Argon 2) System

Now affordable for GI Endoscopy Centers and convenient for physician and patient schedules

SS-200E FEATURES (ESU):

- Fully microprocessed
- Cut and Coagulation Features: **Ecut**, pure cut, blend 1, blend 2, blend 3, spray and bipolar coagulation
- **Ecut** function with 4 homeostasis levels
- Power setting through insulated keys
- Digital display for easy visualization of maximum power levels in each mode
- Digital memory storage capacity for power levels in every mode (pur cut, blend, coagulation and bipolar)
- Automatic voltage selection
- Monitoring system for patient plate circuit
- Permits connection to Canady argon plasma coagulator
- Controlled through manual pencil or footswitch

CANADY PLASMA COAGULATOR ARGON 2 FEATURES:

- Digital display for argon gas flow
- Adjustable argon flow between 0 and 5 l/min
- Adjustment increments for gas flow 0.5 l/min
- Includes trolley cart and accessory kit
- Argon Plasma Cut
- Permits use of footswitch operated Canady pencils and catheters



SS-200E Microprocessed Electrosurgical Unit with Canady Plasma™ Coagulator Argon 2

Total System (SS-200E/Argon 2)

Introductory Pricing:

\$15,995.00

Canady Plasma Probes for Flexible Endoscopy

<u>Quantity</u>	<u>Price</u>
10 (per box)	\$1650.00

Straight-Fire Probes
Side Fire Probes

The single flexible endoscopic Canady Plasma Straight Fire and Side Fire Probes are specifically designed for GI Endoscopy, and Ambulatory Surgery Centers.

To learn more about this system call Parks Medical Today
509-448-9112